

DATA ARTICLE

Open Access



Monitoring of Paraquat in soya products intended for animal feed

Franziska Heydebreck

Abstract

Even though paraquat is considered to be one of the most toxic herbicides, it is one of the most used ones. Paraquat-based preparations are heavily used across many important agricultural countries and for a lot of types of crops such as soya. As soya is used in large quantities in the feed production, paraquat residues in soya products are a potential threat to livestock. The aim of this study was to investigate whether soya products intended for animal feed contain residues of paraquat. Therefore, an optimized Quick Polar Pesticides (QuPPE) Method was developed. In total, 174 samples were analyzed via LC-MS/MS using zwitterionic hydrophilic interaction liquid chromatography (ZIC-HILIC). The results showed that paraquat was detected in 66 % of the samples above 0.008 mg/kg. Residues were mainly present in soya bean meals that were imported from Brazil or the USA and cultivated using genetically modified soya. Around 5 % of the samples were not compliant with the European regulation on maximum residue levels of pesticides in or on food and feed. The results of this study support the view that paraquat residues play a significant role in soya products intended for animal feed and, thus, need to be controlled on a regular basis. However, for enforcement actions reliable processing factors for paraquat in soya products or maximum residue levels for processed soya products are urgently required.

Keywords: Paraquat, Pesticides, Herbicides, Soya, Animal feed, European Union, Brazil, GMO, QuPPE

Introduction

Paraquat is widely considered to be one of the most toxic herbicides currently used. It has been shown that the presence of paraquat can cause harmful effects on terrestrial and aquatic ecosystems (Babatunde et al., 2001; Leboulanger et al., 2009; Muangphra et al., 2014) and human health (Wesseling et al., 2001; Tanner et al., 2011). For example, as paraquat tends to accumulate mainly in the lung, various pulmonary effects, such as wheeze and shortness of breath, have been associated with paraquat exposure of farmers (Hoppin et al., 2002; Schenker et al., 2004). Furthermore, it has been reported that chronic paraquat exposure might contribute to the development of Parkinson's disease in humans (Dinis-Oliveira et al., 2006; Costello et al., 2009; Tanner et al.,

2011). Nevertheless, it is still one of the best-selling herbicides worldwide used for a huge variety of crops, including maize, rice and soya beans (Syngenta Crop Protection AG, 2019).

Paraquat is known for its very effective and fast-acting behavior in killing weeds. It may also be used as a desiccant to hasten maturity of a crop to prepone the harvest date (Boudreaux and Griffin, 2011). However, when paraquat is applied shortly before harvest, residues may be leaved on plant and crop and, thus, may pose a threat to consumers. In the European Union (EU), the authorization of paraquat-based preparations is suspended since 2007. Thus, it is currently not allowed to sell and apply paraquat in the member states of the EU (EuG, 2007). However, paraquat is still authorized and extensively used across many important agricultural countries, including the USA and Brazil (Syngenta Crop Protection AG, 2019).

Correspondence: Franziska.Heydebreck@laves.niedersachsen.de
Feed Institute, Lower Saxony State Office for Consumer Protection and Food Safety, Heckenweg 6, 21680 Stade, Germany



© The Author(s). 2021 **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit <http://creativecommons.org/licenses/by/4.0/>. The Creative Commons Public Domain Dedication waiver (<http://creativecommons.org/publicdomain/zero/1.0/>) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

Chemically seen, paraquat (1,1'-dimethyl-4,4'-bipyridinium) is a cation associated twice with chloride as anion. It belongs to the group of bipyridyl herbicides and quaternary ammonium salts. Paraquat is a very polar substance with a high solubility in water of 620 g/L (National Center for Biotechnology Information, 2020). Because of its ionic nature and high polarity, the analysis of paraquat has been reported in the literature as difficult and challenging. Numerous analytical methods and techniques have been suggested, including gas chromatography (Kawase et al., 1984), capillary electrophoresis (Mallat et al., 2001), immunoassay (Garcia-Febrero et al., 2014), voltammetry (Souza et al., 2006), and spectrophotometry (Shivhare and Gupta, 1991). Liquid chromatography coupled with mass spectrometry via electrospray ionization has become the most common approach, although several problems with the chromatographic (e.g. limited retention and poor peak shape) and mass spectrometric (e.g. various types of quasi-molecular ions and multiple reaction monitoring transitions) behavior of paraquat are described.

The soya bean is the most important oil plant worldwide, with a global oil plant cultivation area of 46 % in 2017 (Ovid, 2018). Soya is a key factor in animal feeding due to its high protein content and profitable cultivation. Soya beans and soya products intended for animal feed are mainly imported from non-EU countries into the EU. In 2017, for example, 14.4 million tons and 20.8 million tons soya beans and soya bean meal, respectively, were imported to the EU (Ovid, 2018). The majority of soya beans was imported from the USA and Brazil (> 70 %), whereas about 87 % of the soya bean meal originated from Brazil and Argentina (Ovid, 2018). In those countries, almost 100 % of the cultivation area is cultivated with genetically modified soya like glyphosate-resistant varieties (Ovid, 2018). However, the intensive use of glyphosate has led to the development of many glyphosate-resistant weeds over the past years (Powles, 2008; Perry et al., 2016; Syngenta Crop Protection AG, 2019). Because paraquat is an important alternative non-selective herbicide, it is very likely that paraquat-based formulations are used as harvest aids to soya beans in those countries. Consequently, imported crops can contain residues of paraquat that might be above permitted tolerances and also be a significant route to a prolonged exposure of both animals and humans to paraquat.

The purpose of this study was to investigate whether soya products intended for animal feed contain residues of paraquat. Furthermore, differences in concentration levels between the various soya products and regions of origin as well as genetically modified soya and non-

genetically modified soya should be evaluated and discussed.

Materials and methods

Samples

Samples were chosen and collected by official feed control inspectors. All sampling sites were located in the German state Lower Saxony and covered a broad spectrum of facilities, including various manufacturers of compound feeding stuffs and the port of Brake on the Weser River, which is an important import cargo spot for feed stock in Europe. The aim was to take samples evenly balanced between various soya products, regions of origin and import countries, and between genetically modified and non-genetically modified material. A particular focus was set on soya bean meal because it is used most for the production of compound feeding stuffs. In total, 174 soya products were collected during the years 2018 and 2019 (Table 1).

Table 1 Overview of the sample material

	No. of samples
Soya product^a	
Soya bean meal (partly steam-heated or dehulled)	144
Soya beans (untreated, steam-heated or toasted)	11
Soya bean protein concentrate	10
Soya bean expeller	7
Soya bean hulls	2
Genetically modified soybeans	
Genetically modified organism (GMO)	94
Non-GMO	68
Unknown	12
Region of origin	
Brazil	88
Argentina	11
Russia	11
India	6
USA	4
Ukraine	2
PR China	2
Austria	2
Canada	1
Lithuania	1
North America	1
Norway	1
Unknown ^b	44

^aDetailed definitions of the various soya products are given in the official catalogue of feed materials (European Commission, 2017a)

^bEspecially for samples from manufacturers of compound feeding stuffs, the region of origin and the point of entry into the EU had to be marked as "Unknown" because there is no obligation to label where feed materials were cultivated

Chemicals

Acetonitrile and methanol (both HPLC gradient grade; $\geq 99.9\%$) were purchased from Honeywell (Seelze, Germany). Formic acid (ACS/Reag. Ph. Eur.; 98–100%), ammonium formate (for LC-MS; $\geq 99.0\%$), and hydrochloric acid fuming 37% (for analysis) were purchased from Merck (Darmstadt, Germany). Ultra-pure water of type 1 grade (American Society for Testing and Materials, ASTM) was supplied by a Sartorius arium system (Göttingen, Germany). Reference standards of paraquat dichloride (CAS-No. 1910-42-5; purity $97.5\% \pm 1.0\%$) and paraquat diiodide-d6 (purity 96.5%) were purchased from LGC Standards (Wesel, Germany). Standard stock solutions of paraquat dichloride and paraquat diiodide-d6 were prepared at approx. 1 mg/mL in methanol and were stored in a freezer at $-20\text{ }^{\circ}\text{C}$. Before use, they were placed in an ultrasonic bath for 10 min to obtain complete dissolution and they were allowed to reach room temperature afterwards. From the standard stock solutions, working standard solutions with concentrations of approx. 10 $\mu\text{g/mL}$, respectively, were prepared and handled according to the standard stock solutions. The working standard solutions were used to prepare fresh daily spiking and calibration solutions with concentrations of approx. 0.4 $\mu\text{g/mL}$ and 0.08 $\mu\text{g/mL}$ in methanol.

Sample preparation

Samples were extracted following the *Quick Polar Pesticides (QuPPE) Method* developed by the EU Reference Laboratory for Single Residue Methods (Anastassiades et al., 2016). However, modifications were made to the QuPPE-Method and are discussed in detail in the [Results and Discussion](#) section.

Samples were quartered and milled to a final particle size of 0.5 mm using an ultra-centrifugal mill. 5 g \pm 0.05 g were weighed into a 50 mL polypropylene tube, spiked with 25 μL of the paraquat diiodide-d6 working standard solution (10 $\mu\text{g/mL}$) and blended. 25 mL of the extraction solvent (methanol/0.5 M hydrochloric acid, 60/40, v/v) were added, blended using a Vortex mixer, and extracted for 15 min using a horizontal shaker. The extract was then centrifuged at 4369 x g and $10\text{ }^{\circ}\text{C}$ for 15 min. Approx. 1 mL of the supernatant was filtrated through a 0.2 μm regenerated cellulose membrane syringe filter (CHROMAFIL[®] Xtra RC-20/13, Macherey-Nagel) and diluted 1:2 with a mixture of methanol and ultra-pure water (60/40, v/v). 5 μL of this solution was injected into the LC-MS/MS system.

Instrumental analysis

Samples were analyzed according to the QuPPE-Method (Anastassiades et al., 2016), but likewise with modifications (cf. section [Results and Discussion](#)). All samples

were analyzed through an LC-MS/MS system using a 1260 Infinity system coupled to a 6470-mass spectrometer (both Agilent Technologies). Analysis was performed using positive electrospray ionization, with unit resolution in Q1 and Q3, and in multiple reaction monitoring (MRM) mode. MRM transitions were 186 > 174 m/z (quantification) and 186 > 77 m/z (confirmation) for paraquat and 192 > 174 m/z for paraquat-d6. High performance liquid chromatography was performed using a SeQuant ZIC[®]-HILIC column (2.1 \times 150 mm, 5 μm , 200 Å ; Merck) combined with an Opti-Lynx Guard column (2.1 \times 15 mm; Optimize Technologies). The mobile phase consisted of (A) ultra-pure water, added with ammonium formate (0.02 mol/L) and adjusted to pH 3 with formic acid, and (B) acetonitrile. The gradient profile was achieved at a flow rate of 200 $\mu\text{L/min}$ and initiated with an equilibration phase of 20% A for 7 min, which was increased to 80% A within 4 min and hold for 10 min. The column was heated constantly at $40\text{ }^{\circ}\text{C}$.

Quality assurance and quality control

Calibration solutions were injected at the start and end of every sample sequence and bracketing calibration was used for quantification. All calibration graphs were linear and non-weighted. The calibration levels usually ranged from 0 ng/mL to approx. 12 ng/mL ($\hat{=}$ 0,060 mg/kg) (4-point calibration). When higher concentrations were quantified, samples were re-analyzed including one higher calibration point. Concentrations were quantified using solvent based calibration curves obtaining the relative response of the target analyte to the amount of the mass-labeled internal standard. Sample concentrations were calculated on a basis of the paraquat dication, so a conversion factor of 0.598 was used to take into consideration that the reference standard was sold as a salt and with remaining water content. The analytical method was fully validated according to the *Analytical quality control and method validation procedures for pesticide residues analysis in food and feed* (European Commission, 2015; European Commission, 2017b). During the method validation, blank samples were spiked with paraquat at concentrations of 38 $\mu\text{g/kg}$, 15.2 $\mu\text{g/kg}$, and 7.6 $\mu\text{g/kg}$ ($n = 5-10$). For these spiked levels, precision values of 6.5%, 4.1%, and 5.1% and trueness values of 100%, 94%, and 89% were yielded, respectively. The limit of quantification (LOQ) was defined as the lowest spiked level that meets the method performance acceptability criteria (i.e. precision $\leq 20\%$ and trueness between 70 and 120%). Accordingly, the LOQ was set to 0.008 mg/kg. As paraquat tends to interact with glass-surfaces, polypropylene tubes and vials were used at all stages of the sample preparation and analysis. Instrumental blanks

could be reduced by using capillaries made of stainless steel instead of polyether ether ketone (PEEK). During a sample sequence, numerous blank injections were made to minimize the carryover of calibration injections. One procedural blank sample was extracted with every sample batch. In general, procedural blank contaminations were low (\bar{O} 0.001 mg/kg) and originated most probably from instrumental contamination. They were not taken into consideration when calculating sample concentrations. The modified method has been successfully tested (z -score = 0.0) in a ring trial between eight laboratories. The ring trial material was used afterwards as reference material in every sample batch to check trueness and precision (precision = 12 %; trueness = 87 %; n = 33). In addition, various spiked matrices (sunflower seed meal, rape seed meal, barley) were analyzed in every sample batch (precision = 15 %; trueness = 89 %; n = 34).

Results and discussion

Modifications to the QuPPE-Method

As mentioned in the section [Materials and Methods](#), modifications were made to the QuPPE-Method because process efficiencies and matrix interferences varied considerably and unsystematically between various soya bean meal samples when using the standard method protocol. A similar observation was also made for the

analysis of paraquat and diquat in cowpea samples (Pizzutti et al., 2016). First of all, water has not been added to the sample homogenate (10 mL of water is typically given to 5 g of sample) because, for some tested soya bean meals, the extraction efficiencies turned out to be up to 55 % lower when water was added. This finding may be in contrast to the assumption, that the addition of water assists the wetting process of dry samples to make pesticides better available for the extraction (Pizzutti et al., 2012). Moreover, it turned out that the soaking time of the water also seems to be not negligible because the process efficiency (defined according to Buhrman et al., 1996 - includes extraction efficiency and matrix interferences) was shown to be up to 30 % lower when soaking time was 60 min instead of no soaking time, concluding that (1) extraction efficiencies were lower, (2) matrix interferences were higher or (3) a combination of both occurred. Following the standard method protocol for extraction (10 mL of 1 % formic acid in methanol as extraction solution; 15 min extraction at 80 °C), process efficiencies were not satisfying and even around 0 % for some test materials. Low recoveries using the standard method protocol were also shown for the analysis of paraquat in other matrices like wheat and chia seeds (Bauer et al., 2018). Optimal extraction conditions were evaluated using various extraction parameters (Table 2).

Table 2 Parameters tested for the extraction of paraquat from soya bean meal. Parameters that were finally used for the extraction of the samples are indicated in bold

Parameter	Variation
Particle size	1 mm; 0.5 mm
Water addition	Yes (10 mL); None
Soaking time of water	0 min; 15 min; 30 min; 45 min; 60 min
Extraction apparatus	Horizontal shaker ; orbital shaker; ultrasonic bath
Extraction solvent	1 % formic acid in MeOH 5 % formic acid in MeOH 5 % formic acid in MeOH/H ₂ O (45/50, v/v) MeOH/0.1 M HCl (50/50, v/v) MeOH/0.5 M HCl (50/50, v/v) MeOH/0.1 M HCl (60/40, v/v) MeOH/0.5 M HCl (60/40, v/v) MeOH/0.75 M HCl (60/40, v/v) MeOH/1 M HCl (60/40, v/v)
Extraction volume	10 mL; 25 mL
Extraction temperature	Room temperature ; 80 °C
Extraction time	15 min ; 30 min; 60 min; 90 min; 120 min
Soaking time of extraction solvent	None ; 15 min; 30 min; 45 min; 60 min
Dilution of extract solution	Undiluted; 1:10; 1:5; 1:2
Clean-up	Solid phase extraction using Strata® WCX; Oasis® HLB; None

There was a slight trend that samples from Argentina and India yielded worse process efficiencies than those from Brazil. Variations in process efficiencies could be attributed to different protein and carbohydrate contents and compositions or different varieties of soya beans and were not investigated in detail.

In addition to the extraction procedure, the chromatographic analysis was also optimized. Paraquat has, in general, a limited retention and poor peak shape on reversed phase columns, which are commonly used in multi-analyte LC-MS/MS analysis, because of its low hydrophobicity and basic properties. The QuPPE-method provides an alternative approach using the Obelisc R™ column (SIELC Technologies) to achieve an optimum retention and peak shape for paraquat. The Obelisc R column consists of a traditional long hydrophobic chain and negatively charged groups so that Paraquat can be retained by a cation-exchange mechanism (SIELC Technologies, 2015). However, bad chromatograms were partly yielded with the Obelisc R

column depending on the sample material (Fig. 1). Most probably, matrix components could not be separated accurately and, thus, interfered with the analyte and suppressed the analyte signal. In conclusion, it was not possible to use the Obelisc R column for the analysis of paraquat in soya bean meal samples. As paraquat is very polar and ionic, another approach for its analysis is hydrophilic interaction liquid chromatography (HILIC) (Hao et al., 2013). With HILIC, paraquat was analyzed in the past in, for example, human urine, cabbage, lake water, palm oil, cereals, and fruits (Whitehead et al., 2010; Robb and Eitzer, 2011; Zou et al., 2014; Francesquett et al., 2019; Halim et al., 2019; Oulkar et al., 2019). For the analysis of paraquat in soya products, there is, to the best of my knowledge, no single residue method published elsewhere. However, based on information and experiences given in the literature, a zwitterionic stationary phase (SeQuant ZIC®-HILIC, Merck) was tested and showed a good retention, electrospray ionization efficiency, and reproducibility (Fig. 1).




Fig. 1 Chromatogram for paraquat (186 > 171 m/z) analyzed in soya bean meal using an Obelisc R™ column (top) and a SeQuant ZIC®-HILIC column (bottom)


Consequently, all samples were analyzed using the SeQuant ZIC®-HILIC column.

Monitoring results

Figure 2 presents an overview of the findings of paraquat residues, highlighting the significant differences between the various regions of origin and soya products and between GMO and Non-GMO material. Paraquat was detected in 66 % of the samples above the LOQ of 0.008 mg/kg. The majority of the samples that contained residues came from Brazil with 69 %, followed by Argentina with 4 %, and the USA and Russia with 3 %, respectively. In Brazil, Argentina, and the USA, more than 94 % of the cultivation area is cultivated with genetically modified soya (Ovid, 2018). Given this situation, it is hardly surprising that paraquat residues were almost entirely found in samples from those countries because it is well-known that paraquat-based formulations can still be effectively used when, for example, glyphosate-resistant weeds have developed (Syngenta Crop Protection AG, 2019). A considerable proportion of 18 % of the found paraquat residues could not be assigned to a

region of origin so that these data is, unfortunately, rather difficult to interpret. The samples from Austria, Canada, India, Lithuania, PR China, and Ukraine were below the LOQ of 0.008 mg/kg. Those samples were mainly declared as Non-GMO material (except the sample from Canada) and, moreover, the use of paraquat has not been authorized in at least Austria and Lithuania as being members of the EU. Furthermore, 60 % of the samples with paraquat residues were cultivated using genetically modified soybeans, whereas 32 % were cultivated using Non-GMO material. These results agree with the findings of the different regions of origin. Also, the Non-GMO material containing paraquat residues came mainly from Brazil (76 %), followed by Russia (11 %), and Norway (3 %). Figure 2 also reveals that of the different soya products, 88 % of paraquat residues were found in soya bean meal and 9 % in soya bean protein concentrate. Soya beans and soya bean hulls had the lowest proportion of paraquat residues at 2 %, respectively, followed by soya bean expeller with 0 %.

In general, the concentrations ranged from 0.008 to 0.070 mg/kg. As shown in Fig. 3, the concentrations for




most of the samples were below 0.04 mg/kg. Only 7 % of the quantified samples were above 0.04 mg/kg and, thus, above the European maximum residue level for soya beans (cf. section *Conclusions*). The highest paraquat concentration of 0.070 mg/kg was found in one sample GMO-soya bean meal from Brazil. In general, the samples with concentrations above 0.04 mg/kg were mainly soya bean meals that were cultivated (1) in Brazil or the USA and (2) using genetically modified soya. A higher concentration of 0.06 mg/kg was also registered for one sample of soya bean hulls.

Conclusions

For all agricultural products for food and animal feed, pesticide maximum residue levels are set by the European Commission and, thus, harmonized throughout the EU via the Regulation (EC) No 396/2005 (European Commission, 2005). A maximum residue level (MRL) is the highest level of a pesticide residue that is legally tolerated in or on food or feed. According to Article 19 of Regulation (EC) No 396/2005, agricultural products that exceed MRLs must not be processed and/or mixed for dilution purposes with a view to placing them on the market as food or feed or feeding them to animals. MRLs apply to fresh products, but also to processed products, considering that the processing may cause a dilution or concentration of the pesticide. When assessing whether a sample contains a residue which exceeds the MRL, a default expanded measurement uncertainty of $\pm 50\%$ has to be applied. For soya beans, an MRL of

0.02 mg/kg is set, which means that results above 0.04 mg/kg are not compliant with the MRL. For processed soya products like soya bean meal, the enforcement of the MRL is limited. Therefore, so-called processing factors must be considered when checking for compliance. For some pesticide-matrix combinations, processing factors are available via databases, for example via the *European database of processing factors for pesticides in food* (Scholz et al., 2018). Those processing factors are not legally binding. Nevertheless, they are a useful tool for the assessment of pesticide residues in food and feed. For paraquat in soya products, reliable processing factors are not available, which complicates considerably the enforcement of the MRL according to Regulation (EC) No 396/2005. Our results clearly show that paraquat residues in soya products might exceed the MRL for soya beans of 0.02 mg/kg or 0.04 mg/kg taking into account the default expanded measurement uncertainty. Around 5 % of our samples were not compliant with the MRL and, thus, required national enforcement actions. However,

The lack of processing factors, especially for soya bean meal as the most important soya product in the feed industry, is problematic and complicates the legal enforcement of the MRL. For the future, it is desirable to establish (1) reliable processing factors for paraquat in soya products or (2) MRLs for processed products to strengthen the national authorities in the EU in controlling and enforcing the Regulation (EC) No 396/2005.



Abbreviations

EU: European Union; GMO: Genetically modified organism; HILIC: Hydrophilic interaction liquid chromatography; LOQ: Limit of quantification; MRL: Maximum residue level; MRM: Multiple reaction monitoring; PEEK: Polyether ether ketone; QuPPE: Quick Polar Pesticides Method; ZIC-HILIC: Zwitterionic hydrophilic interaction liquid chromatography

Acknowledgements

The author thanks Birte Ristau and Theis Oellrich for their assistance in the laboratory, Jens Künemeyer for useful discussions and proofreading, Daniel Heydebreck for helping to create the figures, and the Working Group for Organic Analysis of the Association of German Agricultural Analytic and Research Institutes (VDLUFA) for scientific exchange and participation in the ring trial.

Authors' contributions

The author(s) read and approved the final manuscript.

Funding

As part of the Lower Saxony State Office for Consumer Protection and Food Safety, the Feed Institute is funded by using public money.

Availability of data and materials

The dataset supporting the conclusions of this article is available in the Zenodo repository via <https://doi.org/10.5281/zenodo.3828295>.

Declarations

Competing interests

The authors declare no competing financial interest.

Received: 31 March 2021 Accepted: 26 July 2021

Published online: 12 September 2021

References

- Anastassiades M, Kolberg DI, Benkenstein A et al (2016) Quick Method for the Analysis of numerous Highly Polar Pesticides in Foods of Plant Origin via LC-MS/MS involving Simultaneous Extraction with Methanol (QuPPE-Method). EU Reference Laboratory for Single Residue Method (EURL-SRM), Stuttgart, Germany, Version 9.2
- Babatunde MM, Oladimeji AA, Balogun JK (2001) Acute Toxicity of Gramoxone to *Oreochromis niloticus* (Trewavas) in Nigeria. *Water Air Soil Pollut* 131(1):1–10. doi:<https://doi.org/10.1023/A:1011959514500>
- Bauer A, Luetjohann J, Rohn S et al (2018) Development of an LC-MS/MS Method for Simultaneous Determination of the Quaternary Ammonium Herbicides Paraquat, Diquat, Chlormequat, and Mepiquat in Plant-Derived Commodities. *Food Anal Methods* 11(8):2237–2243. doi:<https://doi.org/10.1007/s12161-018-1201-6>
- Boudreaux JM, Griffin JL (2011) Application Timing of Harvest Aid Herbicides Affects Soybean Harvest and Yield. *Weed Technol* 25(1):38–43. doi:<https://doi.org/10.1614/WT-D-10-00045.1>
- Buhrman DL, Price PI, Rudewicz PJ (1996) Quantitation of SR 27417 in human plasma using electrospray liquid chromatography-tandem mass spectrometry: A study of ion suppression. *J Am Soc Mass Spectrom* 7(11): 1099–1105. doi:[https://doi.org/10.1016/S1044-0305\(96\)00072-4](https://doi.org/10.1016/S1044-0305(96)00072-4)
- Costello S, Cockburn M, Bronstein J et al (2009) Parkinson's Disease and Residential Exposure to Maneb and Paraquat From Agricultural Applications in the Central Valley of California. *Am J Epidemiol* 169(8):919–926. doi:<https://doi.org/10.1093/aje/kwp006>
- Dinis-Oliveira RJ, Remião F, Carmo H et al (2006) Paraquat exposure as an etiological factor of Parkinson's disease. *Neurotoxicology* 27(6):1110–1122. doi:<https://doi.org/10.1016/j.neuro.2006.05.012>
- EuG (2007) Judgment of 11 July 2007, Kingdom of Sweden/Commission of the European Communities, T-229/04, EU:T:2007:217
- European Commission (2015) Analytical quality control and method validation procedures for pesticide residues analysis in food and feed (SANTE/11945/2015)
- European Commission (2017a) Commission Regulation (EU) 2017/1017 of 15 June 2017 amending Regulation (EU) No 68/2013 on the Catalogue of feed materials
- European Commission (2017b) Analytical quality control and method validation procedures for pesticide residues analysis in food and feed (SANTE/11813/2017)
- European Commission (2005) Regulation (EC) No 396/2005 of the European Parliament and of the Council of 23 February 2005 on maximum residue levels of pesticides in or on food and feed of plant and animal origin and amending Council Directive 91/414/EEC
- Francesquett JZ, Rizzetti TM, Cadaval TRS et al (2019) Simultaneous determination of the quaternary ammonium pesticides paraquat, diquat, chlormequat, and mepiquat in barley and wheat using a modified quick polar pesticides method, diluted standard addition calibration and hydrophilic interaction liquid chromatography. *J Chromatogr A* 1592:101–111. doi:<https://doi.org/10.1016/j.chroma.2018.12.060>
- García-Febrero R, Salvador J-P, Sánchez-Baeza F, Marco M-P (2014) Rapid method based on immunoassay for determination of paraquat residues in wheat, barley and potato. *Food Control* 41:193–201. doi:<https://doi.org/10.1016/j.foodcont.2014.01.008>
- Halim N, Kuntom A, Shinde R, Banerjee K (2019) Determination of Paraquat Residues in Palm Oil by High-Performance Liquid Chromatography with UV and Tandem Mass Spectrometry. *Eur J Lipid Sci Technol* 121(8):1900092. doi:<https://doi.org/10.1002/ejlt.201900092>
- Hao C, Zhao X, Morse D et al (2013) Optimized liquid chromatography tandem mass spectrometry approach for the determination of diquat and paraquat herbicides. *J Chromatogr A* 1304:169–176. doi:<https://doi.org/10.1016/j.chroma.2013.07.033>
- Hoppin JA, Umbach DM, London SJ et al (2002) Chemical Predictors of Wheeze among Farmer Pesticide Applicators in the Agricultural Health Study. *Am J Respir Crit Care Med* 165(5):683–689. doi:<https://doi.org/10.1164/ajrccm.165.5.2106074>
- Kawase S, Kanno S, Ukai S (1984) Determination of the herbicides paraquat and diquat in blood and urine by gas chromatography. *J Chromatogr A* 283:231–240. doi:[https://doi.org/10.1016/S0021-9673\(00\)96258-7](https://doi.org/10.1016/S0021-9673(00)96258-7)
- Leboulanger C, Bouvy M, Pagano M et al (2009) Responses of Planktonic Microorganisms from Tropical Reservoirs to Paraquat and Deltamethrin Exposure. *Arch Environ Contam Toxicol* 56(1):39–51. doi:<https://doi.org/10.1007/s00244-008-9164-z>
- Mallat E, Barzen C, Abuknesha R et al (2001) Fast determination of paraquat residues in water by an optical immunosensor and validation using capillary electrophoresis-ultraviolet detection. *Anal Chim Acta* 427(2):165–171. doi:[https://doi.org/10.1016/S0003-2670\(00\)01016-3](https://doi.org/10.1016/S0003-2670(00)01016-3)
- Muangphra P, Kwankua W, Gooneratne R (2014) Genotoxic effects of glyphosate or paraquat on earthworm coelomocytes. *Environ Toxicol* 29(6):612–620. doi:<https://doi.org/10.1002/tox.21787>
- National Center for Biotechnology Information (2020) PubChem Database. Paraquat, CID = 15939. <https://pubchem.ncbi.nlm.nih.gov/compound/Paraquat>. Accessed 4 May 2020
- Oulkar D, Shinde R, Khan Z, Banerjee K (2019) High throughput residue analysis of paraquat and diquat involving hydrophilic interaction liquid chromatographic separation and mass spectrometric determination. *Food Addit Contam Part A* 36(1):120–130. doi:<https://doi.org/10.1080/19440049.2018.1547424>
- Ovid (2018) Verband der ölsaatenverarbeitenden Industrie in Deutschland e. V. <https://www.ovid-verband.de/positionen-und-fakten/ovid-diagramme/>. Accessed 30 Apr 2020
- Perry ED, Ciliberto F, Hennessy DA, Moschini G (2016) Genetically engineered crops and pesticide use in U.S. maize and soybeans. *Sci Adv* 2 (8). doi:<https://doi.org/10.1126/sciadv.1600850>
- Pizzutti IR, Kok A de, Cardoso CD et al (2012) A multi-residue method for pesticides analysis in green coffee beans using gas chromatography–negative chemical ionization mass spectrometry in selective ion monitoring mode. *J Chromatogr A* 1251:16–26. doi:<https://doi.org/10.1016/j.chroma.2012.06.041>
- Pizzutti IR, Vela GME, de Kok A et al (2016) Determination of paraquat and diquat: LC-MS method optimization and validation. *Food Chem* 209:248–255. doi:<https://doi.org/10.1016/j.foodchem.2016.04.069>
- Powles SB (2008) Evolved glyphosate-resistant weeds around the world: lessons to be learnt. *Pest Manag Sci* 64:360–365. doi:<https://doi.org/10.1002/ps.1525>
- Robb CS, Eitzer BD (2011) The Direct Analysis of Diquat and Paraquat in Lake Water Samples by per Aqueous Liquid Chromatography. *LCGC North Am* 29(1):54–59

- Schenker MB, Stoecklin M, Lee K et al (2004) Pulmonary Function and Exercise-associated Changes with Chronic Low-Level Paraquat Exposure. *Am J Respir Crit Care Med* 170(7):773–779. doi:<https://doi.org/10.1164/rccm.200403-266OC>
- Scholz R, van Donkersgoed G, Herrmann M et al (2018) Database of processing techniques and processing factors compatible with the EFSA food classification and description system FoodEx 2 Objective 3: European database of processing factors for pesticides in food. *EFSA Support Publ* 15(11):1510E. doi:<https://doi.org/10.2903/sp.efsa.2018.EN-1510>
- Shivhare P, Gupta VK (1991) Spectrophotometric method for the determination of paraquat in water, grain and plant materials. *Analyst* 116(4):391–393. doi:<https://doi.org/10.1039/AN9911600391>
- SIELC Technologies (2015) Obelisc™ HPLC Columns - Liquid Separation Cell Technology. https://www.sielc.com/wp-content/uploads/2015/11/SIELC_Obelisc_Intro.pdf. Accessed 14 Apr 2020
- Souza D, De, Machado SAS, Pires RC (2006) Multiple square wave voltammetry for analytical determination of paraquat in natural water, food, and beverages using microelectrodes. *Talanta* 69(5):1200–1207. doi:<https://doi.org/10.1016/j.talanta.2005.12.045>
- Syngenta Crop Protection AG (2019) Paraquat information center. <http://www.paraquat.com/en>. Accessed 7 May 2020
- Tanner CM, Kamel F, Ross GW et al (2011) Rotenone, Paraquat, and Parkinson's Disease. *Environ Health Perspect* 119(6):866–872. doi:<https://doi.org/10.1289/ehp.1002839>
- Wesseling C, Van Wendel De Joode B, Ruepert C et al (2001) Paraquat in developing countries. *Int J Occup Environ Health* 7(4):275–286. doi:<https://doi.org/10.1179/oeht.2001.7.4.275>
- Whitehead RD, Montesano MA, Jayatilaka NK et al (2010) Method for measurement of the quaternary amine compounds paraquat and diquat in human urine using high-performance liquid chromatography–tandem mass spectrometry. *J Chromatogr B* 878(27):2548–2553. doi:<https://doi.org/10.1016/j.jchromb.2009.09.029>
- Zou T, He P, Cao J, Li Z (2014) Determination of paraquat in vegetables using HPLC-MS-MS. *J Chromatogr Sci* 53(2):204–209. doi:<https://doi.org/10.1093/chromsci/bmu041>

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

